

Smartline Turn

Technical Specifications

Pressure ulcer grade	● I ● II ● III ● IV
Max. patient weight	230 kg
Item No.	996850/103719




Control unit

Measurements	approx. 29 x 21 x 11 cm
Weight	2.8 kg
Features	Advanced compressor with automatic weight adjustment, offering alternating pressure, static and care modes with 3 cycle options (5, 10, 15 min.), automatic lateral positioning, safety touchpad lock, memory function for quick setup
Alarm functions	Auto-lock feature for enhanced safety, and alarms for both low pressure and power loss

Mattress

Measurements	200 x 88 x 18 cm
Weight	9 kg
Material cells	TPU
Material cover	Bi-elastic PU/PES, breathable, slip retardant
Material base	Nylon with fastening straps
Features	Cell-on-cell-construction with 20 cells and 3 static head cells for stability, non-bending connecting tubes with transportation safety lock, CPR emergency ventilation function for critical scenarios



The SMARTLINE Turn mattress combines advanced pressure care technology with automatic lateral turning, helping to keep patients safe, comfortable, and well supported.

Its innovative air cell system reduces pressure on vulnerable areas, while integrated side bolsters add stability and reassurance. Caregivers benefit from effortless one-touch turning, saving time and reducing physical strain, while patients enjoy uninterrupted comfort and restful sleep. In addition to continuous alternating pressure therapy, the mattress offers customizable care with individually deflatable air cells for targeted offloading, such as dedicated heel relief. This flexibility supports better posture, enhances comfort, and delivers tailored protection for every patient.

- Auto-sensing technology – Automatically adjusts to user weight
- Manual comfort control – Easily set to softer or firmer settings
- 7" cell-on-cell design – For added support and safety during power failure
- Cold device plug – With extra-long 5m mains cable
- Alternating pressure mode – Choose from 3 cycle time options
- Lateral rotation mode – Select from 3 turning intervals
- Static mode – Maintains constant low pressure
- Care mode – For special care scenarios or patient needs
- CPR function – Quick-release for emergency deflation
- 2:1 control (A/B cells) – Enhanced pressure relief cycling
- Integrated cable management – Keeps cords secure and tidy

